

# Large Forced Convection Oven—Single Door

## Features

- Ideal for Large Volume Baking Electronic Parts, Drying, Conditioning and Sterilizing
- inner case is SUS 304 and outside body is powder painting .
- Tempered Safety Glass Window to see inside without door opening
- Install double over temperature protection.
- All models are with fan assisted air circulation, air vent on both side t
- Brake wheels easy to move
- Made In Taiwan



▲ DK-800



▲ DK-601

## Specifications

Model No.	DK-601	DK-800
Temp. Range	Ambient +5°C~210°C ( 250°C )	
Temp. Accuracy	±1°C	
Temp. Control	PID	
Temp. Display	LED digital display	
Temp. Sensor	K-TYPE	
Heater	4000 W	
Volume	270 L	480 L
Timer	99 h / 59 min	
Keep Warm	Rock wool fiber	
Shelf (Adjustable)	4	
Power Voltage	220V 50/60Hz 1PH	
Inside Dimension (WxDxH)	600 x 500 x 900 mm	800 x 600 x 1000 mm
OutsideDimension (WxDxH)	780 x 620 x 1520 mm	980 x 720 x 1620 mm
Weight	130 kg	200 kg

Power Voltage : 110V/220V 50Hz/60Hz