

Forced Convection Oven

Features

- Inside material is SUS304 and outside body is powder painting
- Touch-screen panel
- Transparent tempered glass observation window.
- Motor's axle of stainless steel
- Double Overheat temperature protector
- Using the mineral wool for keep warm
- Built-in circulating fan for temperature stability and with vent can let the steam to evaporate
- Made In Taiwan



◀ DK-600



▲ DK-30D



▲ DK-500

Specifications

Model No.	DK-30D/DK-30DT	DK-500D/DK-500DT	DK-600D/DK-600DT
Temp. Range	Ambient +5°C~210°C		
Temp. Accuracy	±1°C		
Temp. Control	PID		
Temp. Display	LED digital display		
Temp. Sensor	K-TYPE		
Heater	500 W	1200 W	1600 W
Volume	30 L	80 L	150 L
Timer	D series are without Timer / DK series are with timer		
Shelf (Adjustable)	2		
Inside Dimension (WxDxH)	330 x 300 x 310 mm	500 x 400 x 400 mm	600 x 500 x 500 mm
Outside Dimension (WxDxH)	430 x 395 x 570 mm	600 x 510 x 720 mm	700 x 660 x 820 mm
Weight	22 kg	47 kg	73 kg

● We accept "customized order" to high temperature Power Voltage : 110V/220V 50Hz/60Hz