

Hybridization shaking oven

Features

- Inside material is SUS304 and outside body is powder painting
- Touch-screen panel
- Speed and temperature with digital display
- Built-in circulating fan for temperature stability
- Full accessories are available ..
- Orbital(bi-directional) / Rotation / Rocking motion
- Made In Taiwan



▲ OH-800D



▲ OH-801D



◀ OH-AC-H



◀ OH-AC-1.5



◀ OH-AC-15



◀ OH-AC-50



▲ AC-210-96



Specifications

Model No.	OH-800D	OH-801RD	OH-801DS	OH-801RK
Shaking	Rotation+orbital	Rotation	Orbital	Rocking
Temp. Range	Ambient +5°C~80°C			
Temp. Accuracy	±0.1°C			
Temp. Control	PID			
Temp. Display	LED			
Rotor Speed	6~50 r.p.m	6~50 r.p.m	-	-
Rotor Capacity	8 pcs (Max)	Supplied 8 pcs (Max)	-	-
Shaker Speed	10~100 r.p.m	-	20~200 r.p.m	6~60 r.p.m
Shaking Width	25 mm	-	-	-
Shaker Capacity	Carry 6 kg			
Timer	999 h / 999 min / 999 sec			
Shaking Plate	310 x 210 mm		310 x 210 mm	
Inside Dimension (WxDxH)	350 x 250 x 425 mm		350 x 255 x 250 mm	
Outside Dimension (WxDxH)	440 x 375 x 660 mm		440 x 375 x 410 mm	
Weight	33 kg		23 kg	

Power Voltage : 110V/220V 50Hz/60Hz