

Orbital Shaker

Features

- Designed for continuous-operation orbital shaking
- Three-point bearing transmission is adopted suitable for heavy loads and continuous duty.
- Rubber non-slip mat to prevent slipping
- Build-in timer is up to 999/h,m,s.and alarm when time out
- Uses brushless AC motor with spring
- Made in Taiwan



▲ TS-505



▲ TS-500

▲ TS-500A

Specifications

Model No.	TS-500	TS-505	TS-500A
Panel Type	Analogue		
Motion	Orbital		Rocking
Shaking Speed	20~250 r.p.m		6~60 r.p.m
Timer	30 m		
Shaking Width	25 mm		35 mm
Shaking Plate	330 x 280 mm		
Overall Dimension (WxDxH)	330 x 280 x 170 mm	330 x 280 x 370 mm	330 x 280 x 250 mm
Weight	11 kg	12 kg	10 kg

Power Voltage : 110V/220V 50Hz/60Hz