

Orbital shaker

Features

- LED 2-screen digital display actual and setting r.p.m.
- Designed for continuous-operation orbital shaking L
- Differential clamps available
- Three-point bearing transmission is adopted suitable for heavy loads and continuous duty.
- Rubber non-slip mat to prevent slipping
- Build-in timer is up to 999/h,m,s.and alarm when time out
- Uses brushless AC motor with spring
- Made In Taiwan



▲ TS-505D



▲ TS-500D



▲ TS-500AD

Specifications

Model No.	TS-500D	TS-505D	TS-500AD
Panel Type	Digital		
Motion	Orbital		Rocking
Rate of Shaking	20~250 r.p.m	20~150 r.p.m	6~60 r.p.m
Timer	999 h / 999 min / 999 sec		
Shaking Width	25 mm		35 mm
Shaking Plate	330 x 280 mm		
Outside Dimension (WxDxH)	330 x 280 x 170 mm	330 x 280 x 370 mm	330 x 280 x 250 mm
Weight	11 kg	12 kg	10 kg

Power Voltage : 110V/220V 50Hz/60Hz