

Orbital Shaker—Flask Clamp Platform

Features

- Touch-screen panel, LED 2-screen digital display
- Actual and setting speed
- With Digital speed controller and timer
- Shaking smoothly and quietly
- Suitable for carrying heavy load
- Long service time
- Memory of power cut , uses brushless AC motor
- platform can accommodate a mix of flask clamps
- Made In Taiwan



▲ TS-580D

▲ TS-560D

▲ TS-520D



Specifications

Model No.	TS-520D	TS-560D	TS-580D
Control Panel	Digital		
Shaking Speed	20~250 r.p.m		
Shaking Width	25 mm	25 mm	50 mm
Timer	999 h / 999 min / 999 sec		
Shaking Plate	400 x 300 mm	600 x 480 mm	960 x 600 mm
Outside Dimension (WxDxH)	400 x 300 x 220 mm	600 x 480 x 220 mm	965 x 610 x 220 mm
Weight	16 kg	31 kg	98 kg

Capacity for flasks

Flask Clamp	125ml	250ml	500ml	1000ml	2000ml	3000ml	4000ml	5000ml
TS-520D	20	12	6	4	2	1	1	1
TS-560D	48	30	20	12	8	6	6	4
TS-580D	96	60	40	24	15	15	12	8

Power Voltage : 110V/220V 50Hz/60Hz