

# Orbital Shaker—Double Platform

## Features

- Touch-screen panel. LED 2-screen digital display
- Actual and setting speed
- Shaking smoothly and quietly
- Suitable for carrying heavy duty
- With variable speed, timer control and continuous Operation
- Memory of power cut ,uses brushless AC motor
- Differential clamps available
- Double platform for large capacity
- Msde In Taiwan



▲ TS-525D

▲ TS-565D

▲ TS-585D



## Specifications

Model No.	TS-525D	TS-565D	TS-585D
Control Panel		Digital	
Shaking Speed		20~150 r.p.m	
Shaking Width		25 mm	
Timer		999 h / 999 min / 999 sec	
Shaking Plate	400 x 300 2 pcs	580 x 460 2 pcs	960 x 600 2 pcs
Outside Dimension (WxDxH)	400 x 300 x 500 mm	600 x 480 x 500 mm	960 x 600 x 500 mm
Weight	20 kg	60 kg	107 kg

## Capacity for flasks

Flask Clamp	125ml	250ml	500ml
TS-525D	40	24	12
TS-565D	96	60	40
TS-585D	192	120	80

Power Voltage : 110V/220V 50Hz/60Hz