

Refrigerated Circulating Baths

Features

- Circulator to unify temperature inside
- Accurate temperature control with protector
- PID controlled speed and temperature with LED digital display
- Compressor will be delayed to start working if power cut
- Using environmental protection cold-media R134a
- Made with SUS 304 stainless steel body with powder coating
- Pump of pressurize sold separately for your option
- Fitted with device for exterior circulator and drainage stopcock for easy cleaning
- Sensor is using PT-100
- BL-720D series have wheels can move easily



▲ Y-K-BL



▲ BL-720D



▲ BL-710D

Specifications

Model No.	BL-720D		BL-710D	
Model	BL-720D-0	BL-720D-20	BL-710D-0	BL-710D-20
Temp. Range	0°C~100°C	-20°C~100°C	0°C~100°C	-20°C~100°C
Temp. Accuracy	±0.1°C			
Temp. Control	PID			
Temp. Sensor	PT-100 Ω			
Temp. Display	LED			
Volume	20 L		10 L	
Refrigerator	1/2 HP	3/4 HP	1/6 HP	1/4 HP
Circulation Pump	10 L / min			
Heater	1600 W		1200 W	
Bath Opening	300 x 210 mm		200 x 150 mm	
Inside Dimension (WxDxH)	300 x 360 x 180 mm		200 x 300 x 150 mm	
Outside Dimension (WxDxH)	420 x 480 x 800 mm		300 x 400 x 680 mm	
Weight	50 kg		30 kg	

● We accept "customized" order to lower temperature

Power Voltage : 110V/220V 50Hz/60Hz