

## ► Refrigerated Cabinet

Power Voltage : 110V/220V 50Hz/60Hz

### Features

- PID digital controller, touch screen panel, digital temperature setting, automatic defrosting system
- With glass doors, triple-pane design, excellent thermal insulation effect, no water, rendering the best display effect
- Base design of the automatic water evaporation, easy to use
- Base with 360 degree turning wheel and adjustable feet, easy to move, bearing capacity is strong



Specifications	LR-310	LR-320	LR-330
Temp. Range		0°C~20°C	
Temp. Accuracy		±1°C	
Temp. Control		PID	
Temp. Display		LED Display	
Volume	500 L	940 L	1450 L
Shelf (Adjustable)	5	10	15
Door	1	2	3
Refrigerator	1/3 HP	1/2 HP	1 HP
	Using environmental protection cold media R134a		
Safety Device	The compressor will be delayed to start working if power cut and with timer for defrost		
Illumination	40 W x 1	40 W x 1	40 W x 3
Material	Steel plate / Coating		
Inside Dimension (WxDxH)	550 x 640 x 1430 mm	1110 x 650 x 1280 mm	1748 x 650 x 1280 mm
Outside Dimension (WxDxH)	640 x 800 x 2065 mm	1260 x 800 x 2065 mm	1878 x 800 x 2065 mm
Weight (Approx.)	112 kg	165 kg	265 kg



## ► -20C Upright Freezer

### Specifications

SKU	FFU-14M5HW	FFU-21M7HW	FFU-14F5HW	FFU-17F5HW	FFU-21F5HW
Liter	403	583	401	473	575
Lock		*	*	*	*
Light	*	*	*	*	*
Shelves	3	4	3	4	4
Basket	1	1	1	1	1
Remark	Defrost : Manual	Defrost : Manual	Frost Free	Frost free Adjustable Temp. control High-temperature alarm	Frost free Door ajar alarm Control lock
Temp. Range	164 x 82 x 68 cm	178 x 82 x 73 cm	150 x 71 x 73 cm	164 x 82 x 68 cm	178 x 82 x 73 cm