

Cooled Incubator

Features

- Automatic high and low temperature alarm to prevent the temperature from going too high and dropping too low
- Built-in circulating fan for temperature stability
- Humidity System, illumination and Timer can be install if required
- The compressor will be delayed to start working if power cut and with timer for defrost. Using environmental protection cold-media R134a
- Internal Power socket
- Touch-screen panel
- Made In Taiwan

NEW

Glass Window (TSC)



▲ LE-519RD

▲ LE-550RD

▲ LE-400RD

▲ LE-2509RD

▲ LE-509RD

▲ LE-80RD

Specifications

Model No.	LE-519RD	LE-550RD	LE-400RD	LE-509RD	LE-2509RD	LE-80RD
Temp. Range	0°C~70°C					
Temp. Accuracy	±0.1°C					
Temp. Control	PID					
Temp. Display	LED digital display actual and setting temperature					
Temp. Sensor	PT-100 Ω					
Heater	900 W	1200 W	1200 W	900 W	900 W x 2	900 W
Volume	590 L	550 L	400 L	150 L	175 L x 2	80 L
Shelf (Adjustable)	4	4	4	2	2	2
Inside Material	Strengthen white porcelain			SUS 304		
Refrigerator	1/3 HP	1/3 HP	1/3 HP	1/4 HP	1/4 HP	1/4 HP
Inside Dimensin (mm) (WxDxH)	680 x 600 x 1450	680 x 580 x 1400	600 x 500x1350	600 x 500 x 500	700 x 500 x 500	500 x 400 x 400
Outside Dimension (mm) (WxDxH)	815 x 800 x 1940	760 x 820 x 1860	670 x 740 x 1820	680 x 700 x 825	780 x 690 x 1720	560 x 550 x 670
Weight	90 kg	90 kg	90 kg	78 kg	150 kg	50 kg

● Available in custome made designs

Power Voltage : 110V/220V 50Hz/60Hz